


# Quantitative Analysis of Mixtures using Synchrotron XANES Spectroscopy



Canadian Light Source Inc.



# Funding Acknowledgements



**NRC-CMRC**



Saskatchewan  
Industry and  
Resources



Western Economic  
Diversification Canada  
Diversification de l'économie  
de l'Ouest Canada

Alberta

INNOVATION AND SCIENCE

ALBERTA HERITAGE  
FOUNDATION FOR  
MEDICAL RESEARCH



The University of  
Western Ontario



SaskPower



Some of the work presented here was conducted at the NSLS and SRC Use of the National Synchrotron Light Source, Brookhaven National Laboratory, was supported by the U.S. Department of Energy, Office of Science, Office of Basic Energy Sciences, under Contract No. DE-AC02-98CH10886.

This work is based upon research conducted at the Synchrotron Radiation Center, University of Wisconsin-Madison, which is supported by the NSF under Award No. DMR-0084402.

# What is a Synchrotron?

A **synchrotron** is a source of brilliant light that can be used to view the microstructure of materials.




## Linac → Booster ring → Electron Storage Ring

- High voltage DC (MeV) is applied to a tungsten cathode releasing electrons into an ultra high vacuum tube;
- The electrons are fed into a linac where they are accelerated by radiofrequency;
- The linac feeds the electrons into a 103 m booster ring where they receive a boost in energy (250 MeV to 2.5 GeV);
- An injection system transfers the electrons to the 171 m main ring to reach the target circulating current (200 mA).

## What does it do?

- emits highly collimated electromagnetic radiation from subatomic charges orbiting at ultrarelativistic energies in a storage ring with diameters of tens to several hundreds of meters.

# Synchrotron schematic



- high photon flux
- broad spectral range
- high polarization
- natural collimation
- small spot size
- stability
- ring structure allows multiple users

## Electromagnetic Spectrum





# Canadian Light Source Mission:

<http://www.lightsource.ca>

*To advance scientific and industrial capabilities ...*



## Phase 1 beamlines

- XAS (3.5 - 40 keV)
- STXM (205-2000 eV)
- IR (mid, far): (450-6000/10-4000  $\text{cm}^{-1}$ )
- Soft x-rays (SGM/PGM): (5.5-250/200-1900 eV)
- PX1 (protein crystallography): (6.5-18.5 keV)

- only synchrotron in Canada
- 200 million CDN\$ investment
- 2.5 to 2.9 GeV synchrotron operating at 200 mA
- 170.88 m circumference


## main experimental hall




# CLS: 2004/2005

---

- Current Machine parameters
  - Stored electron current: 200 mA
  - Lifetime (1/e) @ 100 mA: 10.5 hours
  - All other design parameters met
- Two beamlines currently being commissioned
- Phase 1 operational by spring 2006



*The Official Opening Ceremony - symphony at the synchrotron Oct. 22, 2004*



*The Prime Minister of Canada May 31, 2004*



*The National News from on top of the Booster ring Oct. 21, 2004*



*The Queen of England visits May 19, 2005*



# Where is Saskatoon, Saskatchewan?



Saskatchewan is characterized by prairie and rolling grassland in the south and lakes and forest in the north.

## Saskatoon


- Founded 1883
- Population - 226 000
- Altitude - 482 m, 52° N
- Climate
  - Precip. rain: 10 in/year
  - snow: 4 in/year
  - Summer: 75 F
  - Winter: -2 F
- Economy





agriculture, mining



# Beamlines



Science is done at beamlines extending up to 50 meters away from the source magnet






# Synchrotrons and Environmental Science



## What x-ray techniques are used for environmental science?






- ❖ X-Ray Fluorescence and imaging XRF (1 ppm)
- ❖ X-ray Absorption Spectroscopy (EXAFS, XANES) XANES (1 ppm)  
EXAFS (10 ppm)
- ❖ x-ray diffraction and scattering (including surfaces)


## What questions are we trying to answer?

- ❖ Will the contaminant bind with soil constituents or is it mobile?
  - ❖ will it get into the groundwater
- ❖ Is the species stable or labile?
- ❖ What is the compounds **toxicity**? **bioavailability**?

# XAS: Origin



# Soils and Sediments



XAS spectroscopy has become a powerful tool in environmental chemistry because of these characteristics,

- ❖ any state of matter, almost any sample - **wet**, dry, **amorphous**, etc.;
- ❖ **element and oxidation state** selective detection;
- ❖ minimal sample preparation
  - ❖ ***in-situ*** studies (no chemical artifacts from extraction);
- ❖ can be **bulk or surface** sensitive
- ❖ **low concentrations** (~ppm and up)

# Metal Interactions at Interfaces




## Adsorption

## Co-precipitation

## Surface precipitation

Mineral  
surface  
layers








EXAFS can distinguish different types of bonding and these suggest different chemical stabilities.

It is also possible to distinguish between different adsorption modes - monodentate, bidentate, etc.



# Mo XANES: a local probe

## Differentiating between Mo oxides




spectra are distinct due to differing amounts of *p/d* orbital mixing

Mo K-edge [20000 eV, [Kr]4d<sup>5</sup>5s<sup>1</sup>]

# Speciation and Mapping of Cr in Hanford Sediments






- 340 million Ci of radionuclides in tank waste.
- 1 million curies have leaked to the soil and groundwater.
- estimated cost of cleanup of the entire US weapons complex is between 189 and 265 billion dollars.





- Cr(III) can be oxidized in the field by Mn solids in soils.
- Cr(VI) can be reduced in the field by soil organic matter, sulfides, ferrous iron - biotite, magnetite, etc.

# Elemental Mapping - Sample SX-108-8A




Sample
41-09-39-6AB
41-09-39-7ABC
SX-108-8A
SX-108-7A

Cr(VI) (%)
58
69
63
27




# XANES

(X-ray Absorption Near-edge Structure)



As K-edge [11867 eV,  $[Ar]3d^{10}4s^24p^3$ ]






# Nickel and Sulfur Speciation in Aerosols

Nickel speciation is important for determining possible toxicity




## Sulfur K-edge XANES

- S K-edge XANES indicates that sulfate is the dominant sulfur species in the air sample.



## Nickel K-edge XANES




- Linear combination fitting of the Ni K-edge XANES indicates a  $\text{NiO}/\text{NiSO}_4$  mixture



## LC fit results from NiO/NiSO<sub>4</sub> mechanical mixtures

conc. (mole %)		LC Fits (mole %)	
NiO	NiSO <sub>4</sub>	NiO	NiSO <sub>4</sub>
10.8	89.2	5.8	94.2
27.8	72.2	18.0	82.0
39.3	60.7	28.3	71.7
51.6	48.4	43.8	56.2
61.0	39.0	50.2	49.7
66.3	33.7	57.6	42.4
76.8	23.2	70.2	29.8
83.1	16.9	82.1	17.9
90.5	9.5	100	0
96.9	3.1	100	0
99.1	0.9	100	0



# LC fitting: XANES





## 3 component mixtures

Actual composition (mole %)			XANES LC Fits (mole %)		
NiO	NiS	NiSO <sub>4</sub> .6H <sub>2</sub> O	NiO	NiS	NiSO <sub>4</sub> .6H <sub>2</sub> O
88.1	8.4	3.5	86.8	9.3	3.9
71.9	21.0	7.0	74.7	17.5	7.8
56.9	26.3	16.7	58.9	23.8	17.3
41.3	47.5	11.2	41.1	47.7	11.2
37.5	29.8	32.7	38.1	27.3	34.6
27.0	19.6	53.3	26.6	20.7	52.7
12.6	84.0	3.4	11.3	85.1	3.6

- LC fits performed on mechanical mixtures;
- Linear regression may offer better results in certain cases but applications are limited and real samples may cause difficulties.



# Linear Regression from specific peak features





# Quantification - Strategies

---



- XANES in energy space, EXAFS in k-space, first-derivative absorption spectra,
- typical approach,
  - principal component analysis (PCA) to determine the number of components in the sample
    - » IND function developed to indicate confidence
  - least squares linear combination fitting
    - » uses spectra of pure reference compounds which are combined to give the compositional fractions of these components.
- pre-edge subtraction, normalization, matrix effects ?
- method to quantify the likely fitting quality
  - Important for K-edge versus L-edge fitting;
- How reliable are these techniques?





Thanks for your Attention