

The Teck Cominco logo features the word "teck" in a bold, dark blue font and "cominco" in a lighter blue font. A vertical blue line is positioned to the left of the text, and a horizontal blue line extends from the bottom of this vertical line across the width of the text.

teckcominco

Using Synchrotron Radiation to Characterize Arsenical Smelting Products

Trevor M. Bergfeldt

Research Chemist, Chemistry & Environment
Teck Cominco Research, Teck Cominco Metals Ltd.
Trail, BC


Outline

- Overview of Teck Cominco & Trail Ops
- Teck Cominco and synchrotron radiation
- Results of an As characterization study
 - #2 Baghouse Dust
 - Fe-As-Sb autoclave residue
- Potential applications
- Summary & Conclusions

Teck Cominco

- Canadian-based natural resource group
- Principal activities: mining, smelting, refining
- 2 Refineries
 - Trail, BC
 - Cajamarquilla, Peru
- Mines in Canada, USA and Peru
(Zn, Pb, Cu, Mo, Au, coal)

Trail Operations

- Largest integrated Zn-Pb smelting & refining complex of its type in the world
 - Capacities: 290,000 t/y Zn, 120,000 t/y Pb
 - Ag, In, GeO_2 , Bi, Cd, Au, low alpha Pb
 - sulphuric acid, ammonium sulphate, SO_2 , elemental S, copper arsenate, copper sulphate, sodium antimonate
 - High level of environmental awareness
- 

Teck Cominco & SR

- Teck Cominco & CLS members of CAMIRO
 - organization that promotes and manages collaborative research within mining industry
- CLS sponsored a demo project at APS for CAMIRO members (March 2002)
 - demonstrate the use of SR-based techniques to mining/metallurgical industry (As and Se)
 - option to attend demo and bring samples

XANES & XAFS

- X-ray absorption near edge structure
 - measures characteristic absorption of an X-ray with an electron in a particular core energy level
 - very sensitive to oxidation state
 - main absorption band for As: $4p \leftarrow 1s$
- X-ray absorption fine structure
 - interaction with neighbouring atoms
 - can “sit” on absorption band of an element and “look out” at local environment

Demonstration at APS

- Advanced Photon Source
 - Chicago, IL
 - PNC-CAT beamline
- Sent arsenic samples for analysis
- Samples compared to arsenic standards
 - AsFeS
 - As₂O₃
 - FeAsO₄·2H₂O

#2 Baghouse Dust

- Product from lead refining
- Very fine dust - 60% Sb & 15% As
- Feed to produce copper arsenate and sodium antimonate
- Historically assumed to be mixture of Sb_2O_3 and As_2O_3
- Actual speciation important for present and future processing options

SEM & XRD

- SEM
 - very fine As & Sb oxides ($< 1\ \mu\text{m}$)
- XRD
 - predominantly senarmontite (Sb_2O_3)
 - essentially “amorphous”
 - likely As_2O_3 (virtually identical to Sb_2O_3)
 - possibly As-O-Sb, but lines also identical

SR Results for #2BHD



- XANES (in a matter of seconds)
 - quantitatively As(III)
- XAFS
 - high degree of fine structure (“very crystalline”)
 - 1st NN is O (i.e. As-O)
 - 2nd NN is NOT As (i.e. not As-O-As)
 - due to high [Sb], As-O-Sb most likely
 - mostly Sb_2O_3 with some SbAsO_3

Fe-As-Sb Residue



- Research project to investigate options to stabilize arsenic and antimony for storage
- #2BHD autoclaved with $\text{Fe}_2(\text{SO}_4)_3$ to precipitate $\text{Fe}_x(\text{AsO}_4)_y$
- Residue: 55% Sb, 14% As, 17% Fe
 - SEM: identified Sb-Fe-As-O & Sb-As-Fe-O
 - XRD: possible As_2O_3
 - “unidentified phase with 3 broad peaks, poor crystallinity and possible hydroscopic nature”

SR Results for Residue

- XANES: 30% As^{3+} , 70% As^{5+}
- XAFS: does not match any standards
 - 1st NN is O (i.e. **As-O**)
 - As(V) as AsO_4^{3-} ■ As (III) as unreacted #2BHD
 - 2nd peak is split, so 2nd NN is ...
 - Fe(?) - but not scorodite or amorphous FeAsO_4
 - possibly AsO_4^{3-} adsorbed on ferrihydrite
 - Sb(?) - matches 2nd peak of #2BHD
 - 30% unreacted starting material (#2BHD)



Fourier Transform of XAFS



Summary

■ #2BHD

- As present is quantitatively trivalent
- As^{3+} - O in #2BDH is NOT As_2O_3
- Sb_2O_3 with AsSbO_3 (As-O-Sb)

■ Fe-As-Sb Residue

- ~ 30% As not completely oxidized
- As^{5+} is AsO_4^{3-} , possibly adsorbed on ferrihydrite
- also As-O-Sb, from unreacted starting material

Other applications ...

- Speciation & association of metals in ...
 - ores & concentrates (PGMs, Ag, Au, In, Ge)
 - mill tailings, effluents & other in-process materials
 - ARD
 - soils & sediments
- Speciation of surface coatings on minerals from various portions of flotation circuits
- Map slices of rock/ore/agglomerate
- Surface analysis to look for corrosion, cracks, layered coatings with high resolution

Conclusions

- SR-based techniques valuable tools for mining & metallurgical applications:
 - determine oxidation states
 - determine nearest neighbour associations
- Allows for characterization of material not possible by traditional methods
 - high resolution
- 2 additional projects ongoing with CLS

Acknowledgements

- Teck Cominco Research
 - GGRichards, KMPenney, RWGreen, SJWalden
- Canadian Light Source
 - JNCutler & NChen
- BC Synchrotron Institute
 - CJones

Questions ?

